Checkme O2

Smart Wrist Pulse Oximeter

User Manual

English | Deutsch | Italiano | Español | Français

www.viatomtech.com

Contents

User Manual.....	English 1-20
Benutzerhandbuch.....	Deutsch 21-40
Manuale d'uso.....	Italiano 41-61
Manual de usuario.....	Español 62-82
Manuel de l'utilisateur.....	Français 83-103

User Manual

Contents

1 Introduction	2
2 Using the Monitor	6
3 PC software	10
4 Maintenance	11
5 Troubleshooting	12
6 Specifications	12
7 Electromagnetic Compatibility	15

1 Introduction

1.1 Intended Use

This product is intended to be used for measuring, displaying and storing of pulse oxygen saturation (SpO_2), pulse rate and motion in home or healthcare facilities environment.


Note:

- The data and results provided by this device are for pre-check screening purpose only and cannot be directly used for diagnostic or treatment.
- The data provided by the APP and PC software is not intended for diagnosis or treatment purpose, always consult your doctor for any health condition.

Warnings and Cautionary Advice

- DO NOT squeeze the sensor or apply excessive force on the sensor & cable.

- Do not use this device during MRI examination.
- Never submerge the device in water or other liquids. Do not clean the device with acetone or other volatile solutions.
- Do not place this device in pressure vessels or gas sterilization device.
- Consult your doctor immediately if you experience symptoms that could indicate acute disease.
- Do not self-diagnose or self-medicate on the basis of this device without consulting your doctor. In particular, do not start taking

any new medication or change the type and/or dosage of any existing medication without prior approval.

- Use only cables, sensors and other accessories specified in this manual.
- Prolonged continuous SpO₂ monitoring may increase the risk of undesirable changes in skin characteristics, such as irritation, reddening, blistering or burns.
- Check the SpO₂ sensor application site every 6-8 hours to determine the positioning of the sensor and the circulation and skin sensitivity of the patient. Patient sensitivity varies depending on medical status or skin condition. For patients with poor peripheral blood circulation or sensitive skin, inspect the sensor site more frequently.
- Functional tester cannot be used to assess the accuracy of a SpO₂ sensor or a device.
- This device is designed to determine the arterial oxygen saturation percentage of functional hemoglobin. Factors that may degrade pulse oximeter performance or affect the accuracy of the measurement include the following:
 - excess ambient light
 - excessive motion
 - electrosurgical interference
 - blood flow restrictors
(arterial catheters, blood pressure cuffs, infusion lines, etc.)
 - moisture in the sensor
 - improperly applied sensor
 - incorrect sensor type
 - poor pulse quality
 - venous pulsations

- anemia or low hemoglobin -concentrations
- cardiogreen and other -intravascular dyes
- carboxyhemoglobin
- methemoglobin
- dysfunctional hemoglobin

1.2 Guide to Symbols

Symbol	Description
	Manufacturer
	Date of manufacture
	Serial number
	Unique device identifier
	Indicates that the product should not be discarded as unsorted waste but must be sent to separate collection facilities for recovery and recycling.
	Refer to instruction manual
	Type BF Applied Part
	No alarm system
	MR unsafe
IP22	Indicates that the product is protected against solid foreign objects of 12,5 mm Ø and greater; and protected against vertically falling water drops when enclosure tilted up to 15°.
	Indicates that caution is necessary when operating the device or control close to where the symbol is

	placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences
	Indicates that the product complies with the EU Medical Device Regulations (Regulation (EU) 2017/745)
	Medical device
	Authorized representative in the European community
	UKCA marking
	UK Responsible Person
	Non-ionizing radiation
	Indicates that the marked item or its material is part of a recovery or recycling process.
	Our products and packaging can be recycled, don't throw them away! Find where to drop them off on the www.quefairedemesdechets.fr site (Only applicable for French market).
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation

1.3 Unpacking

- Main Unit × 1
- Ring Sensor × 1

- Cable × 1
- User Manual × 1

2 Using the Monitor

2.1 Charging

Charge the battery before using.

Use the charge cable to charge the battery of device in the USB Port of the computer or with USB charging adapter.

After being fully charged, the device will power off automatically.

2.2 POWER ON/OFF

POWER ON:

Press the button for 1 second to turn on the device.

POWER OFF:

Automatically power off: The device will turn off automatically in 2 minutes if no measurement, no operation or without App connection.

Manually power off: You can press the button for about 2 seconds to turn off.


2.3 TAKE THE FIRST RECORDING

START. Wear the device and the Ring Sensor, press the button to power on. And keep yourself in the quiet environment.

(For the sake of clarity, we recommend that users wear the Checkme O2 watch on their left wrist and put the ring sensor on the thumb. **If it is too tight, try another finger. Avoid being loose.**)

STOP. After the record, take off the Ring Sensor (and the device), the recording will be saved after the countdown. (If the working time is less than 1 minute, the data will not be saved.)

Note:

- The device can also work as a clock and pedometer without ring sensor, and the steps will restart counting after the ring sensor plugged in.
- Please avoid excessive motion for the sensed finger during recording and avoid any strong ambient light condition.
- We recommend using the classic monitor mode in daily environments and the sleep mode when you want to sleep.

2.4 Stop monitoring & sync data

Take off the sensor, the countdown will begin.

During the countdown, if you wear the sensor again, the record will be resumed.

After the countdown, the data will be ready for sync.

Note: The duration of the recording is 1mins~10 hours. And the device can collect and store maximum 4 recordings, when the 5th recording coming, the first recording will be will covered.

Sync data to the ViHealth App

You can sync the data to your app after the countdown or next time after you turn on the device.

Note: Refer to the operation guide in ViHealth app for more details.

2.4.1 How to find the Serial Number?

- 1) Detach the main unit from wristband.
- 2) Flip the main unit over, the series number is printed on the label of product.

Note: The serial number (SN) is on the back of the device

2.5 How to synchronize the time of the device

The time of the Checkme O2 will be automatically synchronized with the network time on your smart device after connected with the app.

Note: Refer to the operation guide in ViHealth app for more details.

2.6 Operation Mode Navigation

The device has two operation modes: Sleep Mode and Monitor Mode. Turn on the device, during the icon blinking, press the button to switch between two modes.

	Sleep Mode	Monitor Mode
Parameters	SpO ₂ , PR, Motion	SpO ₂ , PR, Steps
Screen	Activated by button	Always on
The Sensor Vibrates	When SpO ₂ is lower than threshold set	When SpO ₂ is lower than threshold set:

		Or when step goal is met
Bluetooth	On when the screen lights up	Always on
Report	Sleep report	Monitor report

2.7 Smart Vibration

The vibrator in the ring sensor will be activated when the SpO₂ fall below the preset value (Threshold). The vibration will stop when the SpO₂ recover to normal level or you can press the button to stop it.

You can customize the vibration to help you sleep better via app or PC software.

2.8 Download App

App name: ViHealth

iOS: App Store

Android: Google Play

Compatibility

The device is compatible with iOS versions 9.0+ and Android versions 5.0+. Refer to the operation guide in app for more details..

2.9 Bluetooth Connection

The device Bluetooth will be enabled automatically only when the **screen lights up**.

Note: when the Checkme O2 is pairing with the smart device, keep the screen lights up.


Note: DO NOT PAIR in the settings of your smart device.

2.10 Display

During measurement, press the button, you can switch different displays.

Time, Remaining battery capacity

Sleep Mode

Monitor Mode

SPO2

Blood oxygen saturation

Pulse rate in Sleep Mode

Pulse rate in Monitor Mode

Wear the sensor

Steps

O2 Score, SpO₂ drop times

3 PC software

3.1 Download the PC Software

PC Software: **O2 Insight**

Download from: getwellue.com/pages/pc-software

Install the software on Windows PC or Mac.

3.2 How to connect the Checkme O2 to Windows

1. Turn on the device, connect the device to PC USB port with the supplied cable of data.
2. Open the PC software, download data from the device.
3. With the PC software, you can customize each report.
4. You can view and print sleep report, which can also be exported as PDF or CSV files.

3.3 How to connect the Checkme O2 to Mac

1. Before the connection, make sure your device with any connection with App. Then turn on device, and turn on the Bluetooth of the Mac.
2. Open the “**O2 Insight Pro**” software, and click the “download” on the screen. Then the data report will sync with your MAC.
3. With the optional PC software, you can customize each report.
4. You can view and print sleep report, which can also be exported as PDF or CSV files.

4 Maintenance

4.1 Cleaning

Use a soft cloth moistened with water or alcohol to clean the device surface.

4.2 Battery

To keep the battery in good condition, charge the battery every 6 months when the device is not in use.

5 Troubleshooting

Problem	Possible Cause	Possible Solution
Device does not turn on or no response.	Battery may be low.	Charge battery and try again.
	Unexpected software condition	Press the button for about 10 seconds to reset
	Device might be damaged.	Please contact your local distributor.
The app cannot find the device.	The Bluetooth of your phone is off.	Turn on the Bluetooth in the phone.
	The device Bluetooth is off in Sleep Mode.	Press the button, the Bluetooth will be turned on when the screen lights up.
The device screen displays "Error 1".	Errors occur during data analysis.	Connect the power supply and plug in the charging cable, press and hold the touch key for 3s to reset the hardware.
The device screen displays "Error 4".	Sensor error.	Please re-plug the sensor.

6 Specifications

Classifications	
EC Directive	MDD, 93/42/EEC

	RED, 2014/53/EU	
	ROHS 2.0, 2011/65/EU	
Degree protection against electrical shock	Type BF	
Environmental		
Item	Operating	Storage
Temperature	5 to 40°C	-25 to 70°C
Relative humidity (noncondensing)	10% to 95%	10% to 95%
Barometric	700 to 1060 hPa	700 to 1060 hPa
Degree of dust & water resistance	IP22	
Physical		
Weight	12 g (main unit)	
Display	OLED	
Wireless	Bluetooth 4.0 BLE	
Vibrator	Probe built-in	
Power Supply		
Charge input:	DC 5V ±10%	
Battery type	Rechargeable lithium-polymer battery	
Charge time	2~3 hours	
Battery life	16 hours	
SpO ₂		
Standards	Meet standards of ISO 80601-2-61	
Measurement accuracy verification: The SpO ₂ accuracy has been verified in human experiments by comparing with arterial blood		

sample reference measured with a CO-oximeter. The pulse rate accuracy has been verified by Emulator. Pulse oximeter measurement are statistically distributed and about two-thirds of the measurements are expected to come within the specified accuracy range compared to CO-oximeter measurements.	
SpO2 range	70% to 100%
SpO2 Accuracy (Arms)	80-100%:±2%, 70-79%:±3%
PR range	30 to 250 bpm
PR accuracy	±2 bpm or ±2%, whichever is greater
Wave length	660-940nm
Output power	Red/Infrared: 3mW max. avg.
Pedometer	
Range	0 to 99999 steps
Storage	
Capacity	4 records,10 hours for each
Sleep analysis	
Record parameters	SpO2, pulse rate, motion
Record interval	4s
Drop analysis	Yes
Mobile APP	
iOS	iOS 9.0 or above, iPhone 4s/iPad 3 or above
Android	Android 5.0 or above, with Bluetooth 4.0 BLE
Bluetooth RF	
Frequency range	2.402 – 2.480 GHz
Max RF power	-10 dBm

Durable period	
Expected service life	3 years

7 Electromagnetic Compatibility

The device meets the requirements of EN 60601-1-2.

Warnings and Cautions

- Using accessories other than those specified in this manual may result in increased electromagnetic emission or decreased electromagnetic immunity of the equipment.
- The device or its components should not be used adjacent to or stacked with other equipment.
- The device needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.
- Other devices may interfere with this device even though they meet the requirements of CISPR.
- When the inputted signal is below the minimum amplitude provided in technical specifications, erroneous measurements could result.
- Portable and mobile communication equipment may affect the performance of this device.
- Other devices that have RF transmitter or source may affect this device (e.g. cell phones, PDAs, and PCs with wireless function).

Table 1

Guidance and manufacturer's declaration-electromagnetic emission		
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Pulse Oximeter should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment-

		guidance
RF emissions CISPR 11	Group 1	The Pulse Oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Pulse Oximeter suitable for use in all establishments, including domestic establishments and those directly network that supplies buildings used for domestic purposes.
Harmonic emissions IEC61000-3-2	N/A	
Voltage fluctuations/flicker emissions IEC61000-3-3	N/A	

Table 2

Guidance and manufacturer's declaration-electromagnetic emission			
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Pulse Oximeter should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge(ESD) IEC61000-4-2	±8 kV contact ±15kV air	±8 kV contact ±15kV air	Floors should be wood, concrete or ceramic tile. if floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast	±2kV for	N/A	N/A

transient/ burst IEC61000-4-4	power Supply lines ± 1 kV for input/output lines		
Surge IEC 61000-4-5	± 1 kV line (s) to line(s) ± 2 kV line(s) to earth	N/A	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle <40% UT (60% dip in UT) for 5 cycles <70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	N/A	N/A
Power frequency (50Hz/60Hz) magnetic field IEC61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

Table 3

Guidance and manufacturer's declaration – electromagnetic immunity			
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of The Pulse Oximeter should assure that it is used in such an electromagnetic environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment –guidance
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of The Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$
Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2 \sqrt{P}$ 80MHz to 800MHz $d = 2.3 \sqrt{P}$ 800MHz to 2.5GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). b Field strengths from fixed RF transmitters, as determined by an electromagnetic site

			survey ,a should be less than the compliance level in each frequency range .b Interference may occur in the vicinity of equipment marked with the following symbol. 
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a: Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which The Pulse Oximeter is used exceeds the applicable RF compliance level above, The Pulse Oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating The Pulse Oximeter. b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.			

Table 4

Recommended separation distances between portable and mobile RF communication the equipment
The Pulse Oximeter is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of The Pulse Oximeter can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Pulse Oximeter as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W(Watts)	Separation distance according to frequency of transmitter M(Meters)		
	150kHz to 80MHz $d=1.2 \sqrt{P}$	80MHz to 800MHz $d=1.2 \sqrt{P}$	80MHz to 2,5GHz $d=2.3 \sqrt{P}$
0,01	N/A	0.12	0.23
0,1	N/A	0.38	0.73
1	N/A	1.2	2.3
10	N/A	3.8	7.3
100	N/A	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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Shenzhen Viatom Technology Co., Ltd.
MedNet EC-REP GmbH
MediMap Ltd